

## Directions for use:

Apply directly to the knees, fetlocks, shins and tendons after periods of training or incurred heat. Spread evenly about 6mm thick.

For a more prolonged effect, wrap with moistened brown paper and bandage.

Repeat as necessary after 24 hours.

To remove, hose off with water.

## Nil Withholding period

Does not contravene the rules of racing or equestrian competition within New Zealand.

© Staysound and the TRM symbol are registered trademarks of TRM-Thoroughbred Remedies Manufacturing Ltd, Ireland

## What makes STAYSOUND® so special?

- STAYSOUND® contains Thyme oil, an essential oil and mother natures natural antiseptic, making STAYSOUND safe to use over minor cuts and abrasions.
- STAYSOUND® is specially formulated for a longer cooling action.
- STAYSOUND® can be used right up to race day, and can be continued for as long as necessary.
- STAYSOUND® is easy to apply with no need to wrap or bandage, and is easy to wash off.
- Studies show that leg clays have the best long term cooling effect when compared to ice packs, ice gel boots, branded ice bandages, cold hosing and a proprietary equine ice refrigeration system.
- Aids in the relief of heat, inflammation and shin-soreness in horses.
- Excellent to use on the lower legs after work as a preventative for shin-soreness and tendon damage.

# STAYSOUND®



## Beat the heat! Stay Cool STAYSOUND




Proudly distributed by Domhealth  
PO Box 204093, Highbrook 2161  
Ph: 09 274 7676 | 0800 500 344  
Fax: 09 274 6164  
Email: [meganr@domhealth.co.nz](mailto:meganr@domhealth.co.nz)



# STAYSOUND®

Internationally acclaimed cold clay dressing which cools, soothes and tightens hot, tired legs.



Staysound® can be used as part of the daily routine, to care for and help maintain trouble free legs.

Staysound® has a highly effective cooling action, which is very beneficial for hot, tired legs, especially during demanding times.

The legs are the most vulnerable structure of the horse and yet they must withstand enormous pressures and forces of impact during training, racing and competition. At the end of training, racing or competition cooling the legs effectively will allow your horse to recover more quickly.

Staysound® tightens, giving support and comfort to fatigued legs and helping to reduce inflammation.

Staysound® has a unique formula containing natural thyme oil, known for its non-irritating, extremely high antiseptic properties and specially formulated to produce a longer cooling action. It can be used safely and confidently right up to race day, over minor cuts and abrasions, and can be continued for as long as necessary.

It is easy to apply with no need to wrap and is easy to wash off.

Ready to use, easy to apply with a safe, rapid action cooling. This product is a must for any competitive horse.



**Pack sizes:** 1.5kg, 5kg, 11.35kg and 20kg resealable tubs.

**Ingredients:** Montmorillonites, Aluminium Silicate, Sodium Borate (Borax), Glycerine, Natural Thyme oil, Purified Water

## What trial studies have to say...

In a independent trial study, TRM's Staysound® was found to reliably reduce the lower leg temperature-the only leg clay in the trial to produce lower leg temperatures in all the trial horses.

The study also concluded that leg clays have the best long term cooling effect when compared to ice packs, ice gel boots, branded ice bandages, cold hosing and a proprietary equine ice refrigeration system.

